



2X Taq FroggaMix

#FBTAQM 1ml

Contents:	2X Taq Mix	1ml
	Nuclease-free water	1ml

Store at -20°C

For research use only

Contains total 2 vials.

Description

2X Taq Mix is a premixed, ready-to-use solution containing Taq DNA Polymerase, dNTPs, Mg²⁺ and Reaction Buffer at optimal concentrations for efficient amplification of DNA templates by PCR, as well as an inert loading dye. To prepare the final PCR, only primers and template DNA need to be added. This pre-mixed formulation saves time and reduces contamination due to the fewer pipetting steps required for PCR set up. The mix retains all features of Taq DNA Polymerase.

Taq DNA Polymerase is a thermostable recombinant DNA polymerase derived from thermophilic bacterium *Thermus aquaticus*. The molecular weight is 94 kDa. Taq DNA Polymerase can amplify DNA targets up to 5 kb (simple template). The elongation velocity is 0.9~1.2kb/min (70~75°C). It has 5' to 3' polymerase activity but lacks 3' to 5' exonuclease activity which results in a 3'-dA overhang PCR product.

Applications

- High-throughput PCR.
- Routine PCR with high reproducibility
- Generation of PCR products for TA cloning

Features

- **Convenient** : Just add primers and template DNA
- **High yields** of PCR products with minimal optimization.
- **High efficiency**: saves your time by simplifying the process
- **Reproducible**: lower contamination risk and pipetting error.

Composition of the 2xTaq Mix

0.25U/ul Taq DNA polymerase, 2X PCR buffer, 0.4mM dNTPs, 3.2 mM MgCl₂, 0.02% bromophenol blue.

Taq mix buffer is a proprietary formulation optimized for robust performance in PCR.

Basic PCR Protocol

All solutions should be thawed on ice and gently mixed

1. Add in a thin walled PCR tube on ice:

For a total 50µl reaction volume

Component of sample	Volume	Final concentration
Taq Mix (2X)	25 µl	1X
Forward Primer	variable	0.1-1 µM
Reverse Primer	variable	0.1-1 µM
Template DNA	variable	10 pg-1 µg
Water, nuclease-free	to 50 µl	–

Recommendations with Template DNA in a 50µl reaction volume:

Human genomic DNA	0.1 µg-1 µg
Plasmid DNA	0.5 ng-5 ng
Phage DNA	0.1 ng-10 ng
<i>E.coli</i> genomic DNA	10 ng-100 ng

2. Gently mix the sample and briefly centrifuge to collect all liquid to the bottom of the tube.

3. Overlay the sample with mineral oil or add an appropriate amount of wax. This step may be omitted if the thermal cycler is equipped with a heated lid.

4. Perform PCR using the following thermal cycling conditions.

Initial Denaturation	94°	3 minutes
25-35 Cycles	94°	30 seconds
	55-68°	30 seconds
	72°	1 minutes
Final Extension	72°	10 minutes

5. Incubate for an additional 10 min at 72°C and maintain the reaction at 4°C. The samples can be stored at -20°C until use.

6. Analyze the amplification products by agarose gel electrophoresis and visualize by staining with EtBr or RedSafe™ (cat# 21141). Use appropriate molecular weight standards.

Notes on cycling conditions

- Recombinant Taq DNA Polymerase is the enzyme of choice for most PCR applications.
- The half-life of enzyme is >40 minutes at 95°C.
- The error rate of Taq DNA Polymerase in PCR is 2.2×10^{-5} errors per nt per cycle; the accuracy, (an inverse of the error rate) an average number of correct nucleotides incorporated before making an error, is 4.5×10^{-4}
- Taq DNA Polymerase accepts modified nucleotides (e.g. biotin-, digoxigenin-, fluorescent-labeled nucleotides) as substrates for the DNA synthesis.
- The number of PCR cycles depends on the amount of template DNA in the reaction mix and on the expected yield of the PCR product. 25-35 cycles are usually sufficient for the majority of PCR reaction. Low amounts of starting template may require 40 cycles.

Guidelines for preventing contamination

During PCR more than 10 million copies of template DNA are generated. Therefore, care must be taken to avoid contamination with other templates and amplicons that may be present in the laboratory environment. General recommendations to lower the risk of contamination are as follows:

- Prepare your DNA sample, set up the PCR mixture, perform thermal cycling and analyze PCR products in separate areas.
- Set up PCR mixtures in a laminar flow cabinet equipped with a UV lamp.
- Wear fresh gloves for DNA purification and reaction set up.
- Use reagent containers dedicated for PCR. Use positive

displacement pipettes, or use pipette tips with aerosol filters to prepare DNA samples and perform PCR set up.

- Always perform “no template control” (NTC) reactions to check for contamination

Quality Control

The absence of endodeoxyribonucleases, exodeoxyribonucleases and ribonucleases is confirmed by appropriate quality tests. Functionally tested in amplification of a single-copy gene from human genomic DNA.

Endodeoxyribonuclease Assay

No detectable conversion of covalently closed circular DNA to a nicked DNA was observed after incubation of 25µl Taq Mix (2X) with 1µg pBR322 DNA in 50µl for 4 hours at 37°C and 70°C.

Exodeoxyribonuclease Assay

No detectable degradation of lambda DNA-HindIII fragments was observed after incubation of 25µl Taq Mix (2X) with 1µg digested DNA in 50µl for 4 hours at 37°C and 70°C.

Ribonuclease Assay

0% of the total radioactivity was released into trichloroacetic acid-soluble fraction after incubation of 25µl Taq Mix (2X) with 1µg E.coli [3H]-RNA (40000cpm/µg) in 50µl for 4 hours at 37°C and 70°C.

www.froggabio.com

Tel (416) 736-8325 Fax (416) 645-3399 Toll Free (877) 318-7277
Office 230 Canarctic Drive, Toronto, Ontario M3J 2X8, Canada
E-mail info@froggabio.com

