

pLUG-Prime[®] TA-Cloning Vector Kit II

For fast ligation using TA-cloning vector

RUO

Research Use Only

REF

11063



DESCRIPTION

- pLUG-Prime TA-Cloning Vector Kit II offers a quick, reliable and efficient method for cloning a variety of DNA sequence
- Possible to separate insert DNA from vector using *Eco* RI restriction enzyme
- More accurate results / Accept a wide range of inserts with different sizes
- Two types of ligation buffers provided for your convenience

INTRODUCTION

TA-cloning technology exploits the terminal transferase activity of some DNA polymerases such as *Taq* DNA polymerase and other non-proofreading DNA polymerase. These enzymes preferentially add a 3'-end A-overhang to PCR products. This allows the direct insertion of such PCR products into the prelinearized cloning vector, which has a T-overhang on each 3'-end. This eliminates the need for restriction enzyme digestion of the vector or insert, primers with built-in restriction sites, or specially designed adapters, resulting in a much more efficient and robust cloning procedure. This technique is especially useful when compatible restriction sites are not available for the subcloning of DNA fragments. pLUG-Prime[®] TA-Cloning Vector Kit II contain several engineered restriction-enzyme recognition sites around the TA-cloning site allowing easy restriction analysis of recombinant plasmids or re-cloning to another vector. Especially, the restriction-enzyme recognition sites of pLUG-Prime[®] TA-Cloning Vector Kit II around the TA-cloning site are more general and simple, which is useful in downstream application such as re-cloning to another vector.

❖ Features of the pLUG-Prime[®] TA-Cloning Vector Kit II

• LacZ alpha sequence

The fragment of lacZ alpha sequence in the pLUG-Prime[®] TA-Cloning Vector Kit II is able to complement beta-galactosidase activity. The lacZ alpha sequence reduces the time to screen for positive clones.

• Multiple cloning region

The multiple cloning region is located around the TA-cloning site in the pLUG-Prime[®] TA-Cloning Vector Kit II. The restriction-enzyme recognition sites of pLUG-Prime[®] TA-Cloning Vector Kit II around the TA-cloning site are located in mirror-repeat pattern.

KIT CONTENTS and STORAGE CONDITION

Components	Concentration	Volumes
TA-Cloning Vector II (20 reactions)	25 ng/μl	40 μl
Control insert DNA	10 ng/μl	10 μl
T4 DNA ligase	2 U/μl	20 μl
10X Ligation Buffer A	-	50 μl
10X Ligation Buffer B	-	50 μl
Forward Primer (M13-F)	10 μM	50 μl
Reverse Primer (M13-R)	10 μM	50 μl
Storage Conditions : -20 °C		

It is recommended for the products to be stored at -20 °C. Always avoid multiple freeze-thaw cycles or exposure to frequent temperature changes. These fluctuations can greatly alter stability of product.

CHARACTERISTICS

- High cloning efficiency
- High percentage of true white colony
- Credible blue/white colony selection
- Rapid procedure (rapid ligation)
- Allowing easy re-cloning to another vector
- Allowing convenient sequencing (M13 F/R priming sites)

Cloning into the pLUG-Prime[®] TA-Cloning Vector Kit II

•Optimizing insert to vector ratio

For the ligation reaction, the optimal molar ratio of insert (i.e., PCR product) to vector has to be optimized. We recommend using the 5-10 fold molar excess of PCR product over TA-cloning vector. In some cases, a lower ratio of PCR product to vector may be sufficient for efficient ligation. The pLUG-Prime[®] TA-Cloning Vector Kit II is about 2.73 Kb and the kit suggests adding 50ng (2 μl) of the vector for the ligation reaction.

• Calculated for 50ng vector using the following equation :

$$\text{ng PCR product required} = \frac{50\text{ng} \times \text{PCR product size (bp)}}{\text{Vector size (bp)}} \times \text{molar ratio}$$

Table 1. Guide for the amount of PCR product to use in the ligation reaction

PCR product size (bp)	5-times molar excess (ng)	10-times molar excess (ng)
500	45.8	91.6
1,000	91.6	183.3
1,500	137.5	274.9
4,500	412.4	824.8

Protocol for ligation using the pLUG-Prime[®] TA-Cloning Vector Kit II

1. Centrifuge pLUG-Prime TA-Cloning Vector Kit II and /or PCR DNA tubes to collect contents at the bottom of the tubes.
2. Vortex the ligation buffer vigorously before use.
3. Set up the following items as described below :

Components	Standard control	Positive control
10X Ligation buffer A	1 μl	1 μl
10X Ligation buffer B	1 μl	1 μl
TA-Cloning Vector II	2 μl	2 μl
PCR product	X μl	-
T4 DNA ligase	1 μl	1 μl
Control DNA	-	3 μl
Deionized Water		Up to 10 μl

* To obtain higher white colonies having insert, we recommend to use the gel-purified PCR product.

4. Mix the reactions by pipetting.

5. Incubate the reactions for 5 to 15 min at 22°C. Alternatively, if the maximum of transformants is required, incubate the reactions overnight at 4°C.


6. Implement transformation with appropriate competent cells.**Note :**

- To facilitate the A-addition of PCR product, add 10 minutes of extension at 72°C to the end of the PCR cycling.
- The use of fresh PCR product is strongly recommended. Purification of PCR products prior to ligation will generally result in higher transformation efficiencies. There are several choices to clean up the PCR product. During gel-purification of PCR-products, avoid long exposure to UV-light. To reduce nicking of the DNA, crystal violet staining is an alternative for band visualization.
- PCR products generated using proofreading DNA polymerases such as *Pfu* DNA polymerase can be used in TA-cloning procedures after the addition of a 3'-end A-overhang.
- We recommend using the 5-10-fold molar excess of PCR product over TA-cloning vector. In some cases, a lower ratio of PCR product to vector may be sufficient for efficient ligation.
- If you see precipitate in the ligation buffer, warm it briefly at 37°C to dissolve the salts. It is important to mix the solutions completely before use to avoid localized concentrations of salts.
- The pLUG-Prime® TA-Cloning Vector Kit II is compatible with a wide range of chemically competent cells which are available from iNtRON. In particular, DH5a competent cells are appropriate to blue/white colony selection and exhibit the high transformation efficiency. Use the 3-10µl aliquot of ligation mixture to transform competent cells.

PROTOCOL for Colony PCR

- Pick an isolated colony with a sterile toothpick, and then, dilute the colony in 20µl of DW. Use the appropriate diluted colony as PCR template.

- Set up program of the thermal cycle



* The setting time of 72°C (for example 30 sec) is different according to length of insert DNA. (In general, DNA polymerase can synthesize 1Kb DNA in 1minute.)

- Check the PCR products by using 1% agarose gel :

※ When you check the band in the gel, consider the length of insert DNA between M13 forward primer and M13 reverse primer in the vector.


MAP & MULTIPLE CLONING SITE

Figure 1 : Map and sequence reference points of the pLUG-Prime® TA-Cloning Vector II

* Before the insert is incorporated into the pLUG-Prime® TA-Cloning Vector II, there is only one *HindIII* site and no *BglII* site. After the incorporation, the T and A nucleotide on the insert will complement with the sequence on the vector and generate these two new sites. This merit of pLUG-Prime® TA-Cloning Vector II makes cloning more economical and convenient.

Multiple Cloning Region	434 to 490
LacZ Operator	531 to 548
LacZ gene	511 to 149
Amp ^r gene	2528 to 1671
T7 promoter	402 to 439
M13 forward primer	359 to 375
M13 reverse primer	507 to 528
β-lactamase coding region	1524 to 2528

301	TACGCCAGCT ATGGGGTGA	GCGCAAAGG CCGCTTTCC	GGATGTGCTG CCTACAGCAC	CAAGGCGATT GTTCCGCTAA	AAGTGGGTA TTCAACCCAT
M13 Forward Primer					
351	ACGCCAGGT TGGGTCCA	TTTCCAGTC AAAGGTCAG	ACGACGTGT TCTGCAACA	AAAACGACGG TTTGTCTGCC	CCAGTAATT GGTCACCTAA
T7 Promoter					
401	CTAATACAG CATTATGCTG	TCACATAGG AGTGATATCC	GCAGCTCGG CGCTGAGCC	TACCCAGGG ATGGGCCCGC	AATTCGAAG TTAAGGTTCC
Insert DNA					
451	T T T T T A A -	-A GATCTGGAT T CTAGACCTA	CCCTCTAGA GGGAGATCT	GTCGACCTG CAGCTGGAGG	AGGCAATGCA TCCGTACGTT
M13 Reverse Primer					
491	GCTTGGCGA CGAACCGCT	ATTCTGTGCA TAAGACCAGT	TAGCTGTTTC ATCGACAAG	CTGTGTGAA GAGCAGCTT	TTGTTATCCG AACAAATAGGC


Figure 2 : Multiple cloning site sequence of the pLUG-Prime® TA-Cloning Vector II

SUGGESTIONS

- Avoid multiple freeze-thaw cycles and exposure to frequent temperature changes by preparing single-use aliquots of Ligation Buffers.
- Pfu* DNA polymerase possesses proofreading activity; it does not have the terminal transferase-like activity like Taq DNA polymerase. Ligation reactions using *Pfu* amplified DNA containing no A-tails will result in no positive colonies.
- Methods for increasing the ligation efficiency :
 - A-tailing :

Purified PCR product	X µl
10X PCR buffer	10 µl
10mM dATP	2 µl
Taq	1µl

 - Add deionized water to a final volume of 100ul.
 - Incubate at 72°C for 1 hrs.
 - Purify the A-tailed DNA and use it in the ligation reaction.
 - If the maximum of transformants is required, incubate the reactions overnight at 4°C.
 - The optimal efficiency can be achieved by using a 1:3 molar ratio of vector DNA to the insert DNA.
 - Use competent cells with higher efficiency such as MacCell™ (>10⁸ cfu/µg DNA) series for transformation.
- Using the colony PCR technique, clones can be screened easily and precisely.

TECHNICAL DATA

- Colonies transformed by using pLUG-Prime® TA-Cloning Vector Kit II. (The colony, which is the result of TA-cloning by using pLUG-Prime® TA-Cloning Vector Kit II, was used as the template DNA for PCR.)
- The gel analysis of the PCR products ligated by pLUG-Prime® TA-Cloning Vector Kit II.

