



**FroggaBio**  
scientific solutions


**FREE SHIPPING** orders over \$250

# Novel Juice

**\$79** Buy 2 Get 1 FREE

As low as \$52.67 each

- Ethidium bromide (EtBr) alternative
- Supplied in 6X DNA Loading Buffer, just add to sample
- Perfect for Blue light or UV light detection
- Contains three tracking dyes for visually tracking the DNA migration during the electrophoresis process



Cat#	Description	Price
LD001-1000	Novel Juice DNA Staining Reagent, 1ml	\$79

Toll Free (877) 318-7277

[www.froggabio.com](http://www.froggabio.com)

Prices are valid until March 31<sup>st</sup>, 2021. Price changes and/or early termination of this promotion may occur at any time.

